

Antibodies Product Datasheet

Product Name:	CD40 ligand Antibody
Catalog Number:	C40437
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CD40 ligand
Concentration:	1 mg/ml
Gene Name:	CD40LG CD40L TNFSF5 TRAP
Protein Name:	CD40 ligand
Other Name:	CD40 ligand (CD40-L) (T-cell antigen Gp39) (TNF-related activation protein) (TRAP) (Tumor necrosis factor ligand superfamily member 5) (CD antigen CD154) [Cleaved into: CD40 ligand, membrane form; CD40 ligand, soluble form]
Human Swiss Prot No.:	P29965
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human CD40 ligand
Storage and Stability:	-20°C/1 year