

Antibodies Product Datasheet

Product Name:	CI153 Antibody
Catalog Number:	C40471
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1: 50-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CI153
Concentration:	1 mg/ml
Gene Name:	C9orf153
Protein Name:	CI153
Other Name:	
Human Swiss Prot No.:	Q5TBE3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CI153 at Human
Storage and Stability:	-20°C/1 year