

Antibodies Product Datasheet

Product Name:	COX6B2 Antibody
Catalog Number:	C40518
Reactivity:	Human, Rat, Mouse
Applications:	IHC, ELISA
Dilutions:	IHC-p 1:50-200, ELISA(peptide)1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human COX6B2
Concentration:	1 mg/ml
Gene Name:	COX6B2
Protein Name:	COX6B2
Other Name:	Cytochrome c oxidase subunit 6B2 (Cancer/testis antigen 59) (CT59) (Cytochrome c oxidase subunit VIb isoform 2) (COX VIb-2) (Cytochrome c oxidase subunit VIb, testis-specific isoform)
Human Swiss Prot No.:	Q6YFQ2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human COX6B2
Storage and Stability:	-20°C/1 year