

Antibodies Product Datasheet

Product Name:	COX7B Antibody
Catalog Number:	C40519
Reactivity:	Human, Rat, Mouse
Applications:	IHC, ELISA
Dilutions:	IHC-p 1:50-200, ELISA(peptide)1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human COX7B
Concentration:	1 mg/ml
Gene Name:	COX7B
Protein Name:	COX7B
Other Name:	Cytochrome c oxidase subunit 7B, mitochondrial (Cytochrome c oxidase polypeptide VIIb)
Human Swiss Prot No.:	P24311
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human COX7B
Storage and Stability:	-20°C/1 year