

Antibodies Product Datasheet

Product Name:	COX82 Antibody
Catalog Number:	C40520
Reactivity:	Human, Rat, Mouse
Applications:	IHC, ELISA
Dilutions:	IHC-p 1:50-200, ELISA(peptide)1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human COX82
Concentration:	1 mg/ml
Gene Name:	COX8A COX8 COX8L
Protein Name:	COX82
Other Name:	Cytochrome c oxidase subunit 8A, mitochondrial (Cytochrome c oxidase polypeptide VIII-liver/heart) (Cytochrome c oxidase subunit 8-2)
Human Swiss Prot No.:	P10176
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human COX82
Storage and Stability:	-20°C/1 year