

Antibodies Product Datasheet

Product Name:	CQ108 Antibody
Catalog Number:	C40530
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1: 50-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CQ108
Concentration:	1 mg/ml
Gene Name:	C17orf108 HSD24
Protein Name:	CQ108
Other Name:	
Human Swiss Prot No.:	A8MSI8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CQ108 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year