

Antibodies Product Datasheet

Product Name:	CTGE5 Antibody
Catalog Number:	C40548
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CTGE5
Concentration:	1 mg/ml
Gene Name:	CTAGE5 MEA11 MEA6 MGEA11 MGEA6
Protein Name:	CTGE5
Other Name:	
Human Swiss Prot No.:	O15320
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CTGE5 at Human/Mouse
Storage and Stability:	-20°C/1 year