

Antibodies Product Datasheet

Product Name:	CX7A1 Antibody
Catalog Number:	C40560
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1: 50-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CX7A1
Concentration:	1 mg/ml
Gene Name:	COX7A1 COX7AH
Protein Name:	CX7A1
Other Name:	
Human Swiss Prot No.:	P24310
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CX7A1 at Human/Mouse
Storage and Stability:	-20°C/1 year