

Antibodies Product Datasheet

Product Name:	CXL10 Antibody
Catalog Number:	C40561
Reactivity:	Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from mouse CXL10
Concentration:	1 mg/ml
Gene Name:	CXL10
Protein Name:	CXL10
Other Name:	C-X-C motif chemokine 10 (10 kDa interferon gamma-induced protein) (Gamma-IP10) (IP-10) (Small-inducible cytokine B10) [Cleaved into: CXCL10(1-73)]
Human Swiss Prot No.:	P17515
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Mouse CXL10
Storage and Stability:	-20°C/1 year