

Antibodies Product Datasheet

Product Name:	CYH2 Antibody
Catalog Number:	C40565
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CYH2
Concentration:	1 mg/ml
Gene Name:	CYTH2 ARNO PSCD2 PSCD2L
Protein Name:	CYH2
Other Name:	
Human Swiss Prot No.:	Q99418
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CYH2 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year