

Antibodies Product Datasheet

Product Name:	DG2L7 Antibody
Catalog Number:	C40596
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DG2L7
Concentration:	1 mg/ml
Gene Name:	DGAT2L7P
Protein Name:	DG2L7
Other Name:	
Human Swiss Prot No.:	Q6IED9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DG2L7 at Human
Storage and Stability:	-20°C/1 year