

Antibodies Product Datasheet

Product Name:	DIAP2 Antibody
Catalog Number:	C40599
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DIAP2
Concentration:	1 mg/ml
Gene Name:	DIAPH2 DIA
Protein Name:	DIAP2
Other Name:	Protein diaphanous homolog 2 (Diaphanous-related formin-2) (DRF2)
Human Swiss Prot No.:	O60879
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human DIAP2
Storage and Stability:	-20°C/1 year