

Antibodies Product Datasheet

Product Name:	DYH2 Antibody
Catalog Number:	C40628
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1: 50-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DYH2
Concentration:	1 mg/ml
Gene Name:	DNAH2 DNAHC2 DNHD3 KIAA1503
Protein Name:	DYH2
Other Name:	
Human Swiss Prot No.:	Q9P225
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DYH2 at Human/Mouse
Storage and Stability:	-20°C/1 year