

Antibodies Product Datasheet

Product Name:	ERAS1 Antibody
Catalog Number:	C40651
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1: 50-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ERAS1
Concentration:	1 mg/ml
Gene Name:	ERICH1-AS1 C8orf68
Protein Name:	ERAS1
Other Name:	
Human Swiss Prot No.:	P0C838
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ERAS1 at Human
Storage and Stability:	-20°C/1 year