

Antibodies Product Datasheet

Product Name:	F19A4 Antibody
Catalog Number:	C40685
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human F19A4
Concentration:	1 mg/ml
Gene Name:	FAM19A4 TAFA4
Protein Name:	F19A4
Other Name:	
Human Swiss Prot No.:	Q96LR4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of F19A4 at Human/Mouse
Storage and Stability:	-20°C/1 year