

Antibodies Product Datasheet

Product Name:	Fetuin-A Antibody
Catalog Number:	C40754
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Fetuin-A
Concentration:	1 mg/ml
Gene Name:	AHSG FETUA PRO2743
Protein Name:	Fetuin-A
Other Name:	Alpha-2-HS-glycoprotein (Alpha-2-Z-globulin) (Ba-alpha-2-glycoprotein) (Fetuin-A) [Cleaved into: Alpha-2-HS-glycoprotein chain A; Alpha-2-HS-glycoprotein chain B]
Human Swiss Prot No.:	P02765
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Fetuin-A
Storage and Stability:	-20°C/1 year