

Antibodies Product Datasheet

Product Name:	FGF7 Antibody
Catalog Number:	C40760
Reactivity:	Human, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FGF7
Concentration:	1 mg/ml
Gene Name:	FGF7 KGF
Protein Name:	FGF7
Other Name:	Fibroblast growth factor 7 (FGF-7) (Heparin-binding growth factor 7) (HBGF-7) (Keratinocyte growth factor)
Human Swiss Prot No.:	P21781
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human,Rat FGF7
Storage and Stability:	-20°C/1 year