

Antibodies Product Datasheet

Product Name:	FKBP1BAntibody
Catalog Number:	C40773
Reactivity:	Human, Rat, Mouse
Applications:	IHC, ELISA
Dilutions:	IHC-p 1:50-200, ELISA(peptide)1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FKBP1B
Concentration:	1 mg/ml
Gene Name:	FKBP1B FKBP12.6 FKBP1L FKBP9 OTK4
Protein Name:	FKBP1B
Other Name:	Peptidyl-prolyl cis-trans isomerase FKBP1B (PPIase FKBP1B) (EC 5.2.1.8) (12.6 kDa FK506-binding protein) (12.6 kDa FKBP) (FKBP-12.6) (FK506-binding protein 1B) (FKBP-1B) (Immunophilin FKBP12.6) (Rotamase) (h-FKBP-12)
Human Swiss Prot No.:	P68106
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human FKBP1B
Storage and Stability:	-20°C/1 year