

Antibodies Product Datasheet

Product Name:	FXYD2 Antibody
Catalog Number:	C40789
Reactivity:	Human, Rat, Mouse
Applications:	IHC, ELISA
Dilutions:	IHC-p 1:50-200, ELISA(peptide)1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FXYD2
Concentration:	1 mg/ml
Gene Name:	FXYD2 ATP1C ATP1G1
Protein Name:	FXYD2
Other Name:	Sodium/potassium-transporting ATPase subunit gamma (Na(+)/K(+) ATPase subunit gamma) (FXYD domain-containing ion transport regulator 2) (Sodium pump gamma chain)
Human Swiss Prot No.:	P54710
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human FXYD2
Storage and Stability:	-20°C/1 year