

## Antibodies Product Datasheet

|                               |                                                                                                            |
|-------------------------------|------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | GNG7 Antibody                                                                                              |
| <b>Catalog Number:</b>        | C40828                                                                                                     |
| <b>Reactivity:</b>            | Human, Mouse, Rat                                                                                          |
| <b>Applications:</b>          | IHC, ELISA                                                                                                 |
| <b>Dilutions:</b>             | IHC-p 1:50-200, ELISA(peptide)1:5000-20000                                                                 |
| <b>Modification:</b>          | Unmodified/Total                                                                                           |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                          |
| <b>Purification:</b>          | The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen. |
| <b>Immunogen:</b>             | Synthesized peptide derived from human GNG7                                                                |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                    |
| <b>Gene Name:</b>             | GNG7 GNGT7                                                                                                 |
| <b>Protein Name:</b>          | GNG7                                                                                                       |
| <b>Other Name:</b>            | Guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-7                                          |
| <b>Human Swiss Prot No.:</b>  | O60262                                                                                                     |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                    |
| <b>Specificity:</b>           | This antibody detects endogenous levels of human GNG7                                                      |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                               |