

Antibodies Product Datasheet

Product Name:	GNG8 Antibody
Catalog Number:	C40829
Reactivity:	Human, Mouse, Rat
Applications:	IHC, ELISA
Dilutions:	IHC-p 1:50-200, ELISA(peptide)1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GNG8
Concentration:	1 mg/ml
Gene Name:	GNG8 GNG9 GNGT9
Protein Name:	GNG8
Other Name:	Guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-8 (Gamma-9)
Human Swiss Prot No.:	Q9UK08
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human GNG8
Storage and Stability:	-20°C/1 year