

Antibodies Product Datasheet

Product Name:	GROβ Antibody
Catalog Number:	C40848
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GROβ
Concentration:	1 mg/ml
Gene Name:	CXCL2 GRO2 GROB MIP2A SCYB2
Protein Name:	GROβ
Other Name:	C-X-C motif chemokine 2 (Growth-regulated protein beta) (Gro-beta) (Macrophage inflammatory protein 2-alpha) (MIP2-alpha) [Cleaved into: GRO-beta(5-73) (GRO-beta-T) (Hematopoietic synergistic factor) (HSF) (SB-251353)]
Human Swiss Prot No.:	P19875
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human GROβ
Storage and Stability:	-20°C/1 year