

Antibodies Product Datasheet

Product Name:	IG2R Antibody
Catalog Number:	C40904
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human IG2R
Concentration:	1 mg/ml
Gene Name:	GIG44 PP9974
Protein Name:	IG2R
Other Name:	
Human Swiss Prot No.:	P09565
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IG2R at Human
Storage and Stability:	-20°C/1 year