

Antibodies Product Datasheet

Product Name:	IGS10 Antibody
Catalog Number:	C40910
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1: 50-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human IGS10
Concentration:	1 mg/ml
Gene Name:	IGSF10 CMF608
Protein Name:	IGS10
Other Name:	
Human Swiss Prot No.:	Q6WRI0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IGS10 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year