

Antibodies Product Datasheet

Product Name:	IL-17 Antibody
Catalog Number:	C40919
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human IL-17
Concentration:	1 mg/ml
Gene Name:	IL17A CTLA8 IL17
Protein Name:	IL-17
Other Name:	Interleukin-17A (IL-17) (IL-17A) (Cytotoxic T-lymphocyte-associated antigen 8) (CTLA-8)
Human Swiss Prot No.:	Q16552
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human IL-17
Storage and Stability:	-20°C/1 year