

Antibodies Product Datasheet

Product Name:	IL-17E Antibody
Catalog Number:	C40922
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human IL-17E
Concentration:	1 mg/ml
Gene Name:	IL25 IL17E UNQ3120/PRO10272
Protein Name:	IL-17E
Other Name:	Interleukin-25 (IL-25) (Interleukin-17E) (IL-17E)
Human Swiss Prot No.:	Q9H293
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human IL-17E
Storage and Stability:	-20°C/1 year