

Antibodies Product Datasheet

Product Name:	IL17F Antibody
Catalog Number:	C40923
Reactivity:	Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from mouse IL17F
Concentration:	1 mg/ml
Gene Name:	IL17F
Protein Name:	IL17F
Other Name:	Interleukin-17F (IL-17F) (Cytokine ML-1) (Interleukin-24) (IL-24)
Human Swiss Prot No.:	Q7TNI7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Mouse IL17F
Storage and Stability:	-20°C/1 year