

Antibodies Product Datasheet

Product Name:	IL1RA Antibody
Catalog Number:	C40927
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human IL1RA
Concentration:	1 mg/ml
Gene Name:	IL1RN IL1F3 IL1RA
Protein Name:	IL1RA
Other Name:	Interleukin-1 receptor antagonist protein (IL-1RN) (IL-1ra) (IRAP) (ICIL-1RA) (IL1 inhibitor) (Anakinra)
Human Swiss Prot No.:	P18510
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human IL1RA
Storage and Stability:	-20°C/1 year