

## Antibodies Product Datasheet

|                               |                                                                                                                                                                |
|-------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | IL-31RB Antibody                                                                                                                                               |
| <b>Catalog Number:</b>        | C40932                                                                                                                                                         |
| <b>Reactivity:</b>            | Human, Rat, Mouse                                                                                                                                              |
| <b>Applications:</b>          | WB, ELISA                                                                                                                                                      |
| <b>Dilutions:</b>             | WB 1:1000-2000 ELISA 1:5000-20000                                                                                                                              |
| <b>Modification:</b>          | Unmodified/Total                                                                                                                                               |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                                              |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.                                                      |
| <b>Immunogen:</b>             | Synthesized peptide derived from human IL-31RB                                                                                                                 |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                        |
| <b>Gene Name:</b>             | OSMR OSMRB                                                                                                                                                     |
| <b>Protein Name:</b>          | IL-31RB                                                                                                                                                        |
| <b>Other Name:</b>            | Oncostatin-M-specific receptor subunit beta (Interleukin-31 receptor subunit beta) (IL-31 receptor subunit beta) (IL-31R subunit beta) (IL-31R-beta) (IL-31RB) |
| <b>Human Swiss Prot No.:</b>  | Q99650                                                                                                                                                         |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                        |
| <b>Specificity:</b>           | This antibody detects endogenous levels of Human IL-31RB                                                                                                       |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                                   |