

Antibodies Product Datasheet

Product Name:	IL33 Antibody
Catalog Number:	C40933
Reactivity:	Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from mouse IL33
Concentration:	1 mg/ml
Gene Name:	IL33
Protein Name:	IL33
Other Name:	Interleukin-33 (IL-33) (Interleukin-1 family member 11) (IL-1F11) (Nuclear factor from high endothelial venules) (NF-HEV)
Human Swiss Prot No.:	Q8BVZ5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Mouse IL33
Storage and Stability:	-20°C/1 year