

Antibodies Product Datasheet

Product Name:	IL5 Antibody
Catalog Number:	C40936
Reactivity:	Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from mouse IL5
Concentration:	1 mg/ml
Gene Name:	IL5
Protein Name:	IL5
Other Name:	Interleukin-5 (IL-5) (B-cell differentiation factor I) (Eosinophil differentiation factor) (T-cell replacing factor) (TRF)
Human Swiss Prot No.:	P04401
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Mouse IL5
Storage and Stability:	-20°C/1 year