

Antibodies Product Datasheet

Product Name:	ITF Polyclonal Antibody
Catalog Number:	C40954
Reactivity:	Human, Rat, Mouse
Applications:	IHC, ELISA
Dilutions:	IHC-p 1:50-200, ELISA(peptide)1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ITF Polyclonal
Concentration:	1 mg/ml
Gene Name:	TFF3 ITF TFI
Protein Name:	ITF Polyclonal
Other Name:	Trefoil factor 3 (Intestinal trefoil factor) (hITF) (Polypeptide P1.B) (hP1.B)
Human Swiss Prot No.:	Q07654
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ITF Polyclonal
Storage and Stability:	-20°C/1 year