

Antibodies Product Datasheet

Product Name:	ITM2C Antibody
Catalog Number:	C40957
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ITM2C
Concentration:	1 mg/ml
Gene Name:	ITM2C BRI3 hucep-14 NPD018 PSEC0047
Protein Name:	ITM2C
Other Name:	
Human Swiss Prot No.:	Q9NQX7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ITM2C at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year