

Antibodies Product Datasheet

Product Name:	K1731 Antibody
Catalog Number:	C40965
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1: 50-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human K1731
Concentration:	1 mg/ml
Gene Name:	KIAA1731
Protein Name:	K1731
Other Name:	
Human Swiss Prot No.:	Q9C0D2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of K1731 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year