

Antibodies Product Datasheet

Product Name:	KIAA0649 Antibody
Catalog Number:	C40986
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KIAA0649
Concentration:	1 mg/ml
Gene Name:	PPP1R26 KIAA0649
Protein Name:	KIAA0649
Other Name:	Protein phosphatase 1 regulatory subunit 26
Human Swiss Prot No.:	Q5T8A7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat KIAA0649
Storage and Stability:	-20°C/1 year