

Antibodies Product Datasheet

Product Name:	KIR3.2 Antibody
Catalog Number:	C40989
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KIR3.2
Concentration:	1 mg/ml
Gene Name:	KCNJ6 GIRK2 KATP2 KCNJ7
Protein Name:	KIR3.2
Other Name:	G protein-activated inward rectifier potassium channel 2 (GIRK-2) (BIR1) (Inward rectifier K(+) channel Kir3.2) (KATP-2) (Potassium channel, inwardly rectifying subfamily J member 6)
Human Swiss Prot No.:	P48051
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human KIR3.2
Storage and Stability:	-20°C/1 year