

Antibodies Product Datasheet

Product Name:	KITH_HHV1C Antibody
Catalog Number:	C40994
Reactivity:	
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from KITH_HHV1C
Concentration:	1 mg/ml
Gene Name:	KITH
Protein Name:	KITH_HHV1C
Other Name:	
Human Swiss Prot No.:	P06478
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KITH_HHV1C
Storage and Stability:	-20°C/1 year