

Antibodies Product Datasheet

Product Name:	KR102 Antibody
Catalog Number:	C41007
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KR102
Concentration:	1 mg/ml
Gene Name:	KRTAP10-2 KAP10.2 KAP18-2 KRTAP10.2 KRTAP18-2 KRTAP18.2
Protein Name:	KR102
Other Name:	
Human Swiss Prot No.:	P60368
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KR102 at Human
Storage and Stability:	-20°C/1 year