

Antibodies Product Datasheet

Product Name:	KR105 Antibody
Catalog Number:	C41008
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KR105
Concentration:	1 mg/ml
Gene Name:	KRTAP10-5 KAP10.5 KAP18-5 KRTAP10.5 KRTAP18-5 KRTAP18.5
Protein Name:	KR105
Other Name:	
Human Swiss Prot No.:	P60370
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KR105 at Human
Storage and Stability:	-20°C/1 year