

Antibodies Product Datasheet

Product Name:	KR106 Antibody
Catalog Number:	C41009
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KR106
Concentration:	1 mg/ml
Gene Name:	KRTAP10-6 KAP10.6 KAP18-6 KRTAP10.6 KRTAP18-6 KRTAP18.6
Protein Name:	KR106
Other Name:	
Human Swiss Prot No.:	P60371
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KR106 at Human
Storage and Stability:	-20°C/1 year