

Antibodies Product Datasheet

Product Name:	KRA13 Antibody
Catalog Number:	C41018
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KRA13
Concentration:	1 mg/ml
Gene Name:	KRTAP1-3 B2B KAP1.2 KAP1.3 KAP1.8 KAP1.9 KRATP1.9 KRTAP1.8
Protein Name:	KRA13
Other Name:	
Human Swiss Prot No.:	Q8IUG1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KRA13 at Human
Storage and Stability:	-20°C/1 year