

Antibodies Product Datasheet

Product Name:	KRA54 Antibody
Catalog Number:	C41021
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KRA54
Concentration:	1 mg/ml
Gene Name:	KRTAP5-4 KAP5.4 KRTAP5.4
Protein Name:	KRA54
Other Name:	
Human Swiss Prot No.:	Q6L8H1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KRA54 at Human/Mouse
Storage and Stability:	-20°C/1 year