

Antibodies Product Datasheet

Product Name:	L37A5 Antibody
Catalog Number:	C41028
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1: 50-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human L37A5
Concentration:	1 mg/ml
Gene Name:	LRRC37A5P C9orf29
Protein Name:	L37A5
Other Name:	
Human Swiss Prot No.:	Q49AS3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of L37A5 at Human
Storage and Stability:	-20°C/1 year