

Antibodies Product Datasheet

Product Name:	LC3 A/B Antibody
Catalog Number:	C41031
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LC3 A/B
Concentration:	1 mg/ml
Gene Name:	LC3 A/B
Protein Name:	LC3 A/B
Other Name:	
Human Swiss Prot No.:	Q9H492/Q9GZQ8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat LC3 A/B
Storage and Stability:	-20°C/1 year