

Antibodies Product Datasheet

Product Name:	LIM and SH3 domain protein 1 Antibody
Catalog Number:	C41043
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human LIM and SH3 domain protein 1. at AA range: 90-130
Concentration:	1 mg/ml
Gene Name:	LASP1 MLN50
Protein Name:	LIM and SH3 domain protein 1
Other Name:	LIM and SH3 domain protein 1 (LASP-1) (Metastatic lymph node gene 50 protein) (MLN 50)
Human Swiss Prot No.:	Q14847
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LIM and SH3 domain protein 1
Storage and Stability:	-20°C/1 year