

Antibodies Product Datasheet

Product Name:	LIS1 Antibody
Catalog Number:	C41047
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LIS1
Concentration:	1 mg/ml
Gene Name:	PAFAH1B1 LIS1 MDCR MDS PAFAHA
Protein Name:	LIS1
Other Name:	
Human Swiss Prot No.:	P43034
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LIS1 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year