

Antibodies Product Datasheet

Product Name:	LMA2L Antibody
Catalog Number:	C41051
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LMA2L
Concentration:	1 mg/ml
Gene Name:	LMAN2L VIPL PSEC0028 UNQ368/PRO704
Protein Name:	LMA2L
Other Name:	
Human Swiss Prot No.:	Q9H0V9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LMA2L at Human/Mouse
Storage and Stability:	-20°C/1 year