

Antibodies Product Datasheet

Product Name:	LRRC3 Antibody
Catalog Number:	C41071
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LRRC3
Concentration:	1 mg/ml
Gene Name:	LRRC3 C21orf102 LRRC3A UNQ9233/PRO31982
Protein Name:	LRRC3
Other Name:	
Human Swiss Prot No.:	Q9BY71
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LRRC3 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year