

Antibodies Product Datasheet

Product Name:	LRRT4 Antibody
Catalog Number:	C41074
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LRRT4
Concentration:	1 mg/ml
Gene Name:	LRRTM4 UNQ3075/PRO9907
Protein Name:	LRRT4
Other Name:	
Human Swiss Prot No.:	Q86VH4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LRRT4 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year