

Antibodies Product Datasheet

Product Name:	LSM2 Antibody
Catalog Number:	C41077
Reactivity:	Human, Mouse, Rat
Applications:	IHC, ELISA
Dilutions:	IHC-p 1:50-200, ELISA(peptide)1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LSM2
Concentration:	1 mg/ml
Gene Name:	LSM2 C6orf28 G7B
Protein Name:	LSM2
Other Name:	U6 snRNA-associated Sm-like protein LSm2 (Protein G7b) (Small nuclear ribonuclear protein D homolog) (snRNP core Sm-like protein Sm-x5)
Human Swiss Prot No.:	Q9Y333
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human LSM2
Storage and Stability:	-20°C/1 year